


How to visualize various biodata via Vaa3D



Xin Chen




Institute for Brain and Intelligence

10/18/2021







Biodata – various, complicated, containing valuable information

Various species



Various resolution



Various imaging modalities

Various data types, formats, etc.

Vaa3D – 3D Visualization-Assisted Analysis

Visualization

Various-scale (small and large, e.g. terabyte voxel) multidimensional (3D/4D/5D) bioimage



- Tri-view
- 3D viewer – 3D/4D/5D bioimage
- TeraFly – bioimage with unlimited size
- TeraVR – immersive

Data Types


- 3D color image stacks
 - Tiff stack (.tif, .tiff)
 - Zeiss LSM (.lsm)
 - MRC (used for electron microscopy images) (.mrc)
 - Vaa3D's raw file (.v3draw, .raw, *.v3dpbd, *.vaa3dpbd) (The PBD files are compressed raw files)
 - Any other image formats (e.g. jpeg, PNP, BMP, ...) supported by LOCI Bioinformatics Java library (through a Vaa3D-bioformats plugin, see the plugin pages)
- 5D time series of color image stacks
 - Each time point saved as a separate file (end with suffix like 000.tif, 001.tif, ...)
 - Each time point saved as a single slice of a 3D image stack of whatever formats Vaa3D supports (e.g. tiff, or Vaa3D's raw)
- 3D irregular shaped surfaces
 - Wavefront .OBJ files
 - Vaa3D's surface format (.v3ds)
- 3D neuron structures or any relational data that can be described as a graph
 - SWC file (.swc)
 - Enhanced SWC file (*.eswc)
- 3D point cloud
 - .apo file (a simple CSV format with fixed number of columns)
- 3D landmarks
 - .marker (indeed a simple CSV format)
 - .csv
- Linker files
 - linker files (.ano) that enables opening many data files of different types easily
 - atlas files (.atlas) for managing colocalized/registered image files

Data Types


Vaa3D: main interface




Tiff stack (.tif, .tiff)




SWC file (.swc)



Vaa3D's raw file (.v3draw)




Vaa3D's raw file (*.v3dpbd)



Tri-view

- 3 orthogonal views
- Focus coordinates
- Zoom
- Looking Glass
- Intensity
- Landmarks
- Atlas
- Color channels

.....



3D Viewer

MIP (Maximal intensity projection)

mIP (minimal intensity projection)

Alpha blending

Cross-sectional view

Threshold

Z-thick

Contrast

Channel

Surf/Object

Volume Colormap


Object Manager

Volume Cut

Surface Cut

Rotation, Zoom, Shift



.....




Vaa3D-TeraFly: Overview

TeraFly extends the Vaa3D software to cope with (potentially) **unlimited** sized bioimage even on laptops with a limited amount of system memory (≤ 4 GB) and video card memory (≤ 1 GB)

- fast rendering/visualization of **3/4/5D terabyte-scale** microscopy images
- **instant** zoom-in/out with mouse-scroll
- visualization-assisted **annotation** of 3D objects at different scales
- image format conversion tool (**TeraConverter**) included




Vaa3D-TeraFly: Architecture




Bria, A. et al. Nature Methods, 2016.

Vaa3D-TeraFly: Schema of 3D-5D Formats




Vaa3D-TeraFly: User Interface



TeraConverter and TeraFly

Generate TeraFly data format:
“Advanced”, “Big-Image_Data”,
“TeraConverter”, “Step 1: Import
volume from:”, “Step 2: Convert
volume to:”, “Start”

**Visualize the TeraFly data
format generated:** “Advanced”,
“Big-Image_Data”, “TeraFly”,
“File”, “Open TeraFly Image (3-
5D)”, “RES(.....)”, “OK”, zoom
in/out, rotate
Adjust the control panel such as
shift, resolution, Z-thick, contrast,
Volume Cut, Surface Cut



TeraVR

Integrate immersive and collaborative 3-D visualization, interaction, and hierarchical streaming of teravoxel – scale images.




Wang, Y. et al. Nature Communications, 2019.

TeraVR

TeraVR (single user): see in VR

VR-Farm (multiple users): collaborate in VR

MR-Farm (multiple users)



TeraVR



Summary

Vaa3D – 3D Visualization-Assisted Analysis

Visualize various multidimensional (3D/4D/5D) biodata

- Tri-view
- 3D viewer – 3D/4D/5D bioimage
- TeraFly – bioimage with unlimited size
- TeraVR – immersive

Extended readings

More details can be found in the below papers:

- Peng, H., Ruan, Z., Long, F., Simpson, J.H. & Myers, E.W. V3D enables real-time 3D visualization and quantitative analysis of large-scale biological image data sets. *Nat. Biotechnol.* 28, 348-353 (2010). DOI:10.1038/nbt.1612
- Peng, H., Bria, A., Zhou Z., Iannello, G. & Long, F. Extensible visualization and analysis for multidimensional images using Vaa3D. *Nat. Protoc.* 9, 193-208 (2014). DOI:10.1038/nprot.2014.011
- Peng, H., Tang, J., Xiao, H., Bria, A., et al. Virtual finger boosts three-dimensional imaging and microsurgery as well as terabyte volume image visualization and analysis. *Nat. Commun.* 5, 4342 (2014). DOI:10.1038/ncomms5342
- Bria, A., Iannello, G., Onofri, L. & Peng, H. TeraFly: real-time three-dimensional visualization and annotation of terabytes of multidimensional volumetric images. *Nat. Methods* 13, 192-194 (2016). DOI:10.1038/nmeth.3767
- Wang, Y. et al. TeraVR empowers precise reconstruction of complete 3-D neuronal morphology in the whole brain. *Nat. Commun.* 10, 3474 (2019). DOI:10.1038/s41467-019-11443-y

Thank you!